

# Protocol for RT-PCR

## Isolation of total RNA

- 1) Aspirate the growth media. Add 1ml of Isogen to each dish. Scrape the cell and transfer Isogen to a new 1.5ml tube. Incubate for 5min at RT.
- 2) Add 200 $\mu$ l of chloroform  $\rightarrow$  Shake vigorously for 15sec.  $\rightarrow$  Store for 2~3min at RT.
- 3) Centrifuge at 15000rpm for 15min at 4C.
- 4) Transfer aqueous phase to a new 1.5ml tube. Add 500 $\mu$ l of isopropanol. Store for 5~10min at RT.
- 5) Centrifuge at 15000rpm for 10min at 4C. Remove the supernatant.
- 6) Add 1ml of 70% EtOH. Vortex briefly.
- 7) Centrifuge at 15000rpm for 5min at 4C. Remove the supernatant.
- 8) Dry the pellet briefly and resolve the pellet with 30 $\mu$ l of DEPC H<sub>2</sub>O.
- 9) Measure the concentration of RNA using Spectrophotometer. Make 0.25 $\mu$ g/ $\mu$ l solution with DEPC H<sub>2</sub>O. (\* DEPC : Diethyl polycarbonate)
- 10) Store at -20C

## Reverse transcription and PCR

### A) Digestion of DNA

- 1) RNA (0.25 $\mu$ g/ $\mu$ l in TE)            8 $\mu$ l(=2 $\mu$ g)  
    10x DNase I Reaction Buffer    1 $\mu$ l  
    DNase I                            1 $\mu$ l  
    incubate for 15min at RT.

- 2) Add 1µl of 25mM EDTA and incubate for 10min at 65C to inactivate DNase.

B) Reverse Transcription

**Denature**

- 1) Add the following reagents and vortex.

Oligo(dT) primer	1µl
10mM dNTP	1µl

- 2) Incubate for 5min at 65C. → Put on ice for 1min.

**Anneal**

- 3) Add cDNA synthesis mixture and vortex.

[cDNA synthesis mixture]	
5x First Strand Buffer	4µl
0.1M DTT	1µl
RNase OUT	1µl
RT(Superscript III reverse transcriptase	1µl

**DNA synthesis**

**Terminate reaction**

- 4) 50 C 60 min → 70 C 15 min
- 5) Add 1µl of RNaseH → 37C 20 min (to **remove RNA.**)

c) RT-PCR

- 1) Mix 3µl of ssDNA and 12µl of D.W. (5x dilution).

- 2) Make the PCR reaction mixture in the PCR tube

1

template	1
10x buffer	2.5
MgCl <sub>2</sub> (25mM)	2
dNTPs	1
Taq polymerase	0.2
Forward Primer (50µM)	0.2 µl
Reverse Primer (50µM)	0.2 µl
<u>D.W.</u>	<u>17.9 µl</u>

Total

25  $\mu$ l

7) Set the PCR program

Stage 1: resp=1, 94 C 5min

Stage 2: resp=25,

step1 (**Denature**) 94C 30sec,

step2(**Anneal**) 55C 30sec ,

Step3 (**Extension**) 72C 30sec

Stage 3 : resp=1, 72 C 5 min

stage 4 : reso=1 , 20 C 1 $\infty$

8) Start the reaction

9) Prepare TAE gel

( TAE buffer ( 40mM Tris-acetate , 1mM EDTA , pH 8.3) + 1.5% Agarose)

Apply samples to the gel

Electrophoresis at 100 V for 15 – 20 min.