

Protocol for Western Blotting

Materials

Materials used on Day 3

•Protease inhibitor stock:

1 $\mu\text{g}/\mu\text{L}$ pepstatin A in DMSO

200 μM leupeptin in OG Buffer

200 mM PMSF: Freshly made. Ex) 34.8 mg PMSF in 1 mL isopropanol or DMSO.

•OG buffer: 150 mM NaCl, 20 mM HEPES, 1mM EDTA, pH 7.4

•OG buffer with protease inhibitors (OGB/PI): OG buffer containing 0.5% (v/v) of protease inhibitor stock

*For example, add 5 μL of protease inhibitor stock to 1 mL of OG buffer.

•Lysis buffer (LB): 0.5% Nonidet P-40, 200 mM NaCl, 20 mM HEPES, 1 mM EDTA, pH 7.4

•Lysis buffer with protease inhibitors (LB/PI): lysis buffer containing 0.5% (v/v) of protease inhibitor stock

*For example, add 5 μL of protease inhibitor stock to 1 mL of lysis buffer.

•Agarose-conjugated anti-cMyc antibody: Santa Cruz, Catalogue No. sc-40AC.

•30% Acrylamide/Bis solution, 37.5:1 (2.6% C)

•Running gel buffer: 1.5 M Tris-HCl, pH 8.8

Materials used on Day 5

- **6x sample buffer** (= 6x sample buffer with reducing reagents):
12% SDS, 40% glycerol, 30% β-mercaptoethanol, 300 mM DTT, 120 mM EDTA, 1 mg/mL bromophenol blue, 375 mM Tris-HCl, pH 6.8.
- **Pre-stained standard**: Invitrogen Catalogue No. LC5925
- **10x running buffer**: 250 mM Tris-HCl, 1920 mM glycine, 1% SDS
- **10x transfer buffer**: 250 mM Tris-HCl, 1920 mM glycine
- **10x PBS** (0.1 M sodium phosphate): 13.8 g $\text{NaH}_2\text{PO}_4 \cdot 4\text{H}_2\text{O}$ in 1 L solution, pH adjusted to 7.2 using NaOH.
- **Stacking gel buffer**: 0.5 M Tris-HCl, pH 6.8
- **immuno-blot PVDF Membrane**. Bio-Rad Catalogue 162-0176,

Materials used on Day 6

TBS: 20 mM Tris-HCl, 500 mM NaCl, pH 7.5

T-TBS: 0.2% Tween-20 in TBS

Blocking buffer: PBS with 10% glycerol, 5% milk, 0.2% Tween-20

Antibody: horseradish peroxidase conjugate anti c-Myc mouse monoclonal IgG (Santa Cruz Biochemistry, sc-40 HRP)

Chemiluminescent HRP Substrate: MILLIPORE , Immobilon Western.

Film: Hyperfilm ECL 18x24, Amersham Biosciences, Product Code RPN2103K

Procedure

Day 3 Cell Lysis and immunoprecipitation

(1) Prepare the buffers

Add 0.5% (v/v) of protease inhibitor stock to 0G buffer or lysis buffer to make 0GB/PI or LB/PI, respectively.

Volume of buffer (mL)

0G Buffer: 3.0 ml / dish

Lysis Buffer: 0.75 ml / dish ,

(2) Preparation of lysate

- 1) Wash 293 cell surfaces twice with 1 mL PBS and aspirate.
- 2) Pour 1 mL of LB/PI onto the cells and incubate for 5min on the ice.
- 3) Transfer LB/PI to a micro-centrifuge tube.
- 4) Centrifuge the tube at 15,000 rpm for 10 minutes at 4C and transfer the supernatant to a new tube.

(3) Wash agarose-conjugated anti-cMyc antibody (AC) for use in immunoprecipitation

- 1) Pour 30 μ L of AC into a new tube.
- 2) Add 0.75 mL of LB/PI.
- 3) Spin for a short time at 15,000 rpm, discard the supernatant, and add 0.75 mL of LB/PI.
- 4) Spin for a short time at 15,000 rpm and discard the supernatant.

(4) Immunoprecipitation

Mix cell lysate (Step (2)-2) and washed AC (Step (3)-4) overnight using a rotator at 4C

Day 4 SDS-PAGE and Western Blotting

(1) Preparation of running gel (7.5% acrylamide)

Number of gels	1	2
30% Acrylamide/Bis solution (mL)	3	6
6Running gel buffer (mL)	3	6
50% glycerol	2.4	4.8
dH ₂ O	3.5	7.0
10% SDS	120 μ	240 μ
10% APS	50 μ	100 μ
TEMED	10 μ	20 μ

Total

12 ml

24 ml

- Pour into gel maker.
- Pour 80% Etoh over the surface of the gel and leave for 30min at RT.

(2) Prepare Stacking gel

Number of gels	1	2
30% Acrylamide/Bis solution (mL)	0.65	1.3
Stacking gel buffer (mL)	1.25	2.5
dH ₂ O(mL)	3.05	6.1
10% SDS(μ l)	50	100
10% APS(μ l)	50	100
TEMED(μ l)	5	10

total

5.1

10.2

- Drain the water from the top of the running gel.
- Pour the stacking gel mixture on top of the running gel.
- Put the comb in place.

(3) Prepare reagents

1) Prepare LB/PI with NaCl

1-1) LB + NaCl: Add 58 mg of NaCl (powder) per 1 mL of LB

1-2) Add 5 μ L of protease-inhibitor (PI) stock per 1 mL of the following buffer.

Volume of buffer (mL)

Lysis Buffer (LB): 3.0 ml / sample

LB + NaCl: 1.5 ml / sample

2) Prepare 1x sample buffer

Dilute sample buffer stock using H₂O.

3) Prepare 1x running buffer

Number of tanks	1	2
10x running buffer (mL)	50	100
H ₂ O (mL)	450	900
Total (mL)	500	1000

4) Prepare transfer buffer

Number of membranes	1	2
10x transfer buffer (mL)	100	200
H ₂ O (mL)	700	1400
Ethanol (mL)	200	400
Total (mL)	1000	2000

(Add Ethanol after mixing 10x transfer buffer and H₂O, then fill to 1000 or 2000 mL using H₂O.)

Cool the buffer at 4°C.

5) Prepare blocking buffer

Prepare 100 - 150 mL of blocking buffer for 1 membrane

For 100 mL of buffer:

10x PBS: 10 mL + Glycerol: 10 mL + NaCl: 0.9 g
+ H₂O: 80 mL

→ Mix → Add 5 g of skim milk

→ Mix → Add 0.2 mL of Tween 20. → Fill to 100 mL using DDW.

(4) Prepare the sample

- 1) Spin the cell lysate and agarose mixture for a short time and discard the supernatant.
- 2) Wash the AC with **LB/PI** twice: Add 0.75 mL of LB/PI, mix using the vortex mixer, spin for a short time, and discard the supernatant. Repeat.
- 3) Wash the AC with **LB/PI containing NaCl** twice. Completely aspirate the supernatant.
- 4) Add 40 µL of 1x sample buffer and mix using the vortex mixer.
- 5) Incubate at 95C for 5 minutes.
- 6) Spin for a short time and transfer the supernatant to a new tube.


(5) SDS-PAGE

- 1) Set the gel with the running buffer and apply the marker and the samples.
- 2) Electrophoresis: at 200 V (constant voltage), for approximately 50 min.

(6) Western Blotting

- 1) Cut PVDF membrane into 7.5 x 10 cm², soak in methanol, and then soak in transfer buffer.
- 2) Disassemble the gel, cut the corner of the gel to mark the orientation, and soak the gel in transfer buffer for 5 minutes.

3) Assemble the gel and the membrane.



4) Blot overnight with 30 V (constant voltage) at 4C.
(or 2 hours with 80V)

Day 6 Blocking, reaction with antibody, and detection

- 1) Soak the membrane in TBS for 5 minutes.
- 2) Soak the membrane in blocking buffer for 2 hours.
- 3) React with antibody (for one or two membranes: blocking buffer, 10 mL; Santa Cruz anti cMyc-HRP mouse monoclonal IgG, 10 μ L) for 2 hours.
- 4) Prepare T-TBS (0.2% Tween-20 in TBS):
For one membrane,
TBS: 50 mL
Tween-20: 100 μ L
Wash the membrane with blocking buffer 5 times for 5 minutes each.
- 5) Wash with T-TBS and TBS.
Soak in T-TBS for 10 minutes. → Soak in TBS for 5 minutes. → Soak in T-TBS for 10 minutes. → Soak in TBS for 5 minutes.
- 6) Mix ECL-plus A solution (5 mL) and B solution (125 μ L), and pour over the membrane. Leave at r.t. for 5 minutes.
- 7) Drain ECL-plus, wrap the membrane, and expose the film.
- 8) Develop the film.

Reference

Nakamura K et al. Mol Endocrinol 18: 1461-1470, 2004.