

## Protocol

### <Core curriculum>

#### A. DNA

##### Small scale preparation of plasmid DNA (mini-prep by PureLink™)

(Preparing Cell Lysate)

1. Transfer 1.5 ml of an overnight culture (E.Coli in LB medium) to a microcentrifuge tube.
2. Centrifuge for 1 min at 15000 rpm at 4 degree. Decant the supernatant.
3. Resuspend the pellet with 250 µl of Resuspension Buffer (R3) with RNase A by vortex.
4. Add 250 µl of Lysis buffer (L7) to cells. Mix gently by inverting the capped tube 5 times. **Do not vortex.**
5. Incubate the tube for 5 min at RT. **Do not exceed 5 min.**
6. Add 350 µl of Precipitation Buffer (N4). Mix immediately by gently inverting the capped tube 5 times.  
**Do not vortex.**
7. Centrifuge at 15000 rpm for 10 min to clarify the lysate.

(Binding DNA)

1. Place a Mini column inside a 2-ml Wash Tube.
2. Load the supernatant from above-mentioned step 7 onto the spin column.
3. Centrifuge the column at 15000 rpm for 1 min. Discard the flow-through and place the column back in the tube.
4. Add 500 µl Wash Buffer (W10) with ethanol to the column. Incubate for 1 min at RT.  
Centrifuge at 15000 rpm for 1 min. Discard the flow-through and place the column back in the tube.
5. Add 700 µl Wash Buffer (W9) with ethanol to the column.
6. Centrifuge at 15000 rpm for 1 min. Discard the flow-through.
7. Centrifuge at 15000 rpm for 1 min to remove any residual Wash buffer (W9).  
Discard Wash tube with the flow-through.

(Eluting DNA)

1. Place the spin column in a clean 1.5-ml Recovery Tube.  
Add 75 µl TE to the center of the column.
2. Incubate the column for 1 min at RT.
3. Centrifuge at 15000 rpm for 2 min.
4. The elution tube contains your plasmid DNA. Discard the column.

## **Digestion of plasmid DNA**

Mix reagents described below in a microcentrifuge tube. Incubate it at 37 degree for 60 min.

T7blue-MUC5AC	5 $\mu$ l
10x buffer (M)	3 $\mu$ l
SacI	1 $\mu$ l
SpeI	1 $\mu$ l
H2O	20 $\mu$ l
	<hr/>
	30 $\mu$ l

## **Electrophoresis in an agarose gel**

1. Digested DNA 9  $\mu$ l + 10x loading dye 1  $\mu$ l
2. Apply the samples to an agarose gel.
3. Run a gel at 100V for 30 min.
4. Visualize the DNA under UV.

## **B. Protein**

### **Protein extraction and SDS-PAGE**

1. Cultured cells are lysed and scraped on ice with 300 $\mu$ l lysis buffer
2. Collect the sample and mix vigorously for 20 sec
3. Centrifuge 15000rpm for 10min
4. Collect supernatant and add 100  $\mu$ l 4x sample buffer
5. Boil at 100 degree for 2min
6. Apply the 20  $\mu$ l samples onto an 5-20% acrylamide gel.
7. Run the gel until the dye reaches the bottom.

## **C. Cell Culture**

### **Reviving frozen cells**

1. Place the frozen vial containing NCI-H292 cell at 37C degree waterbath just long enough to thaw contents.
2. Slowly remove the cell suspension and place in 15 ml centrifuge tube that contain 5 ml of pre-warmed RPMI + 10% FBS.
3. Centrifuge for 5 min at 1500 rpm.
4. Aspirate the medium and leave the cell pellet.
5. Gently resuspend the cell pellet with 5 ml of RPMI + 10% FBS.
6. Count the cells in hemocytometer.
7. Transfer the cells to 60 mm dish. Place the dish in a CO<sub>2</sub> incubator.