

# Gene Transfection & Western Blotting Protocol

## Schedule

Day 1 Transfection

Day 2 Changing of medium

Day 3 Cell lysis  
Immunoprecipitation  
Preparation of SDS polyacrylamide gel

Day 4 Electrophoresis  
Blotting

Day 5 Blocking  
Reaction with Antibody  
Washing and Detection

## Materials

### Materials used on Day 1

Calcium Phosphate Transfection Kit: Invitrogen Catalogue No. K2780-01

Falcon tube (Catalogue No. 35-2054: 12 x 75 mm, 5-mL tube)

Cells: 293 cells (human embryonic kidney cells), near confluent in  $\phi$ 60-mm dish

Plasmid: cMyc-tagged human LH receptor/pcDNA3.1, 1  $\mu$ g/ $\mu$ L, purified using Qiagen Plasmid Maxi Kit

### Materials used on Day 3

Protease inhibitor stock:

1  $\mu$ g/ $\mu$ L pepstatin A in DMSO

200  $\mu$ M leupeptin in OG Buffer

200 mM PMSF: Freshly made. Ex) 34.8 mg PMSF in 1 mL isopropanol or DMSO.

OG buffer: 150 mM NaCl, 20 mM HEPES, 1mM EDTA, pH 7.4

OG buffer with protease inhibitors (OGB/PI): OG buffer containing 0.5% (v/v) of protease inhibitor stock

For example, add 5  $\mu$ L of protease inhibitor stock to 1 mL of OG buffer.

Lysis buffer (LB): 0.5% Nonidet P-40, 200 mM NaCl, 20 mM HEPES, 1 mM EDTA, pH 7.4

Lysis buffer with protease inhibitors (LB/PI): lysis buffer containing 0.5% (v/v) of protease inhibitor stock

For example, add 5  $\mu$ L of protease inhibitor stock to 1 mL of lysis buffer.

Agarose-conjugated anti-cMyc antibody: Santa Cruz, Catalogue No. sc-40AC.

30% Acrylamide/Bis solution, 37.5:1 (2.6% C): Bio-Rad, Catalogue No. 161-0158

Running gel buffer: 1.5 M Tris-HCl, pH 8.8

#### Materials used on Day 4

6x sample buffer (= 6x sample buffer with reducing reagents):

12% SDS, 40% glycerol, 30%  $\beta$ -mercaptoethanol, 300 mM DTT, 120 mM EDTA, 1 mg/mL bromophenol blue, 375 mM Tris-HCl, pH 6.8.

Pre-stained standard: Invitrogen Catalogue No. LC5925

10x running buffer: 250 mM Tris-HCl, 1920 mM glycine, 1% SDS

10x transfer buffer: 250 mM Tris-HCl, 1920 mM glycine

10x PBS (0.1 M sodium phosphate): 13.8 g  $\text{NaH}_2\text{PO}_4 \cdot 4\text{H}_2\text{O}$  in 1 L solution, pH adjusted to 7.2 using NaOH.

]

Stacking gel buffer: 0.5 M Tris-HCl, pH 6.8

Bio-Rad Catalogue 162-0176, immuno-blot PVDF Membrane

#### Materials used on Day 5

TBS: 20 mM Tris-HCl, 500 mM NaCl, pH 7.5

T-TBS: 0.2% Tween-20 in TBS

Blocking buffer: PBS with 10% glycerol, 5% milk, 0.2%, 0.2% Tween-20

Antibody: horseradish peroxidase conjugate anti c-Myc mouse monoclonal IgG (Santa Cruz Biochemistry, sc-40 HRP)

ECL Plus: ECL Plus Western Blotting Detection Reagents, Amersham Biosciences, Product Code RPN2132

Film: Hyperfilm ECL 18x24, Amersham Biosciences, Product Code RPN2103K

# Gene Transfection and Western Blotting Protocol

## Day 1 Transfection (Calcium-phosphate-mediated transfection)

### ↓ Prepare solution A and B

Solution A: 150  $\mu$ L

10  $\mu$ g plasmid: x  $\mu$ L (1.0  $\mu$ g/ $\mu$ L plasmid  $\rightarrow$  x = 10)

sterile H<sub>2</sub>O: 150 - 18 - x  $\mu$ L (x = 10  $\rightarrow$  122  $\mu$ L)

2 M CaCl<sub>2</sub>: 18  $\mu$ L (0.24 M CaCl<sub>2</sub> in Solution A)

Solution B: 150  $\mu$ L

2x HBS (Hepes Buffered Saline): 150  $\mu$ L

### ↓ Mix Solution A and Solution B

Using a Pasteur pipette, slowly add solution A dropwise to solution B while bubbling air through solution B with another pipette. Continue until solution A is depleted. This is a slow process that should be performed over 1 – 2 minutes.

### Autopipetter

Pasteur pipette (Provides air to make bubbles in Solution B.)

Pipetteman pipette (Used to slowly add Solution A.)

Solution B in Falcon tube

### ↓ Incubate at room temperature (R.T.) for 15 ~ 30 minutes.

### ↓ Remove the spent medium of the 293 cells in the $\phi$ 60-mm dish and add 5 mL of fresh medium.

### ↓ Add the mixture of Solution A and B dropwise to the media in the dish.

### ↓ Incubate overnight at 37C in a humidified CO<sub>2</sub> incubator.

## Day 2 Changing of medium

Gently change all of the medium.

## Day 3 Cell Lysis

(1) Cell harvest and immunoprecipitation

(2) Preparation of running gel

## (1) Cell harvest and immunoprecipitation

### Step 1: Prepare the buffers

Add 0.5% (v/v) of protease inhibitor stock to 0G buffer or lysis buffer to make 0GB/PI or LB/PI, respectively.

Volume of buffer (mL)

0G Buffer:  $3.0 \times a$  (mL)

Lysis Buffer:  $0.75 \times a + 1.5 \times b$

a: number of dishes

b: number of immunoprecipitation samples

### Step 2: Harvest the cells

- 1) Wash 293 cell surfaces with 0.75 mL of 0GB/PI.
- 2) Pour 0.75 mL of 0GB/PI onto the cells and harvest the cells. Transfer the cell suspension to a micro-centrifuge tube.
- 3) Repeat 2). again.
- 4) Centrifuge the tube at 15,000 rpm for 10 minutes at 4C.

### Step 3: Prepare cell lysate

- 1) Discard the supernatant produced during steps 2-4, add 0.75 mL of LB/PI to the sediment, and place on ice for 30 minutes. During the incubation for 30 minutes on ice, the cells and LB/PI should be mixed using the vortex mixer for a few seconds at 5, 10, 15 and 30 minutes after the start of incubation.
- 2) Centrifuge the tube at 15,000 rpm for 10 minutes at 4C and transfer the supernatant to a new tube.

### Step 4: Wash agarose-conjugated anti-cMyc antibody (AC) for use in immunoprecipitation

- 1) Pour 30  $\mu$ L of AC into a new tube.
- 2) Add 0.75 mL of LB/PI.
- 3) Spin for a short time at 15,000 rpm, discard the supernatant, and add 0.75 mL of LB/PI.
- 4) Spin for a short time at 15,000 rpm and discard the supernatant.

### Step 5: Immunoprecipitation

Mix cell lysate (Step 3-2) and washed AC (Step 4-4) overnight using a rotator at 4C

(2) Preparation of running gel (7.5% acrylamide)

Number of gels	1	2
30% Acrylamide/Bis solution (mL)	2.5	5
Running gel buffer (mL)	2.5	5
10% SDS (mL)	0.1	0.2
H <sub>2</sub> O (mL)	4.9	9.8
Total (mL)	10	20

↓ Mix (Do not make bubbles).

↓ Add the following per 10-mL of gel mixture.

10% APS, 50 $\mu$ L  
TEMED, 5 $\mu$ l

↓ Pour into gel maker.

↓ Pour water over the surface of the gel and leave overnight at RT.

Day 4 SDS-PAGE and Western Blotting

(1) Prepare reagents

1) Prepare LB/PI with NaCl

1-1) LB + NaCl: Add 58 mg of NaCl (powder) per 1 mL of LB

1-2) Add 5  $\mu$ L of protease-inhibitor (PI) stock per 1 mL of the following buffer.

Volume of buffer (mL)

Lysis Buffer (LB): 3.0 x b

LB + NaCl: 1.5 x b

b: number of immunoprecipitation samples

2) Prepare 1x sample buffer

Dilute 6x sample buffer to 1x using H<sub>2</sub>O.

Required amount of 1x sample buffer ( $\mu$ L): 40 x c

c: number of samples

### 3) Prepare 1x running buffer

Number of tanks	1	2
10x running buffer (mL)	50	100
H <sub>2</sub> O (mL)	450	900
Total (mL)	500	1000

### 4) Prepare transfer buffer

Number of membranes	1	2
10x transfer buffer (mL)	100	200
H <sub>2</sub> O (mL)	700	1400
Methanol (mL)	200	400
Total (mL)	1000	2000

(Add methanol after mixing 10x transfer buffer and H<sub>2</sub>O, then fill to 1000 or 2000 mL using H<sub>2</sub>O.)

Cool the buffer at 4°C.

### 5) Prepare blocking buffer

Prepare 100 - 150 mL of blocking buffer for 1 membrane

For 100 mL of buffer:

10x PBS: 10 mL

Glycerol: 10 mL

NaCl: 0.9 g

H<sub>2</sub>O: 80 mL

↓ Mix

↓ Add 5 g of skim milk

↓ Mix

↓ Add 0.2 mL of Tween 20.

↓ Fill to 100 mL using DDW.

### (2) Prepare acrylamide gel

Stacking gel

Number of gels	1	2
30% Acrylamide/Bis solution (mL)	0.33	0.65
Stacking gel buffer (mL)	0.63	1.25
10% SDS (mL)	0.05	0.1
H <sub>2</sub> O (mL)	1.53	3.05
Total (mL)	2.53	5.05



- ↓ Mix (Do not make bubbles).
- ↓ Add the following per 2.53 mL of gel mixture.
  - 10% APS, 50  $\mu$ L
  - TEMED, 1.25  $\mu$ L
- ↓ Drain the water from the top of the running gel.
- ↓ Pour the stacking gel mixture on top of the running gel.
- ↓ Put the comb in place.

### (3) Prepare the sample

- 1) Spin the cell lysate and agarose mixture for a short time and discard the supernatant.
- 2) Wash the AC with LB/PI twice: Add 0.75 mL of LB/PI, mix using the vortex mixer, spin for a short time, and discard the supernatant. Repeat.
- 3) Wash the AC with LB/PI containing NaCl twice. Completely aspirate the supernatant.
- 4) Add 40  $\mu$ L of 1x sample buffer and mix using the vortex mixer.
- 5) Incubate at r.t. for 1 hour and at 95 $^{\circ}$ C for 5 minutes.
- 6) Spin for a short time and transfer the supernatant to a new tube.

### (4) SDS-PAGE

- 1) Set the gel with the running buffer and apply the marker and the samples.
- 2) Electrophoresis: at 200 V (constant voltage), for approximately 1 hour.

### (5) Western Blotting

- 1) Cut PVDF membrane into 7.5 x 10 cm<sup>2</sup>, soak in methanol, and then soak in transfer buffer.
- 2) Disassemble the gel, cut the corner of the gel to mark the orientation, and soak the gel in transfer buffer for 5 minutes.
- 3) Assemble the gel and the membrane.

—

gel holder cassette  
 fiber pad  
 filter paper  
 gel  
 PVDF membrane  
 filter paper  
 fiber pad  
 gel holder cassette

+

4) Blot overnight with 30 V (constant voltage) at 4C.

Day 5 Blocking, reaction with antibody, and detection

- 1) Soak the membrane in TBS for 5 minutes.
- 2) Soak the membrane in blocking buffer for 2 hours.
- 3) React with antibody (for one membrane: blocking buffer, 5 mL; Santa Cruz anti cMyc-HRP mouse monoclonal IgG, 10  $\mu$ L) for 2 hours.
- 4) Prepare T-TBS (0.2% Tween-20 in TBS):  
    For one membrane,  
        TBS: 50 mL  
        Tween-20: 100 $\mu$ L  
    Wash the membrane with blocking buffer 5 times for 5 minutes each.
- 5) Wash with T-TBS and TBS.  
    Soak in T-TBS for 10 minutes. → Soak in TBS for 5 minutes. → Soak in T-TBS for 10 minutes. → Soak in TBS for 5 minutes.
- 6) Mix ECL-plus A solution (5 mL) and B solution (125  $\mu$ L), and pour over the membrane. Leave at r.t. for 5 minutes.
- 7) Drain ECL-plus, wrap the membrane, and expose the film.
- 8) Develop the film.

Reference

Nakamura K et al. Mol Endocrinol 18: 1461-1470, 2004.